Zirlux[®] ST1

DISCO / BLOCCHETTI IN ZIRCONIA

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Zirlux[®] ST1

ZIRKONIUM DISCS / BLOCKS

- Isostatically pressed, pre-sintered

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- isostatisch geperst, pre-gesinterd

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
DE Healthcare Ltd
Unit 9, Kingsthorpe Business Centre,
Studland Road, Kingsthorpe,
Northampton NN2 6NE

Distributed by/Distribuido por/Distribué par/Vertrieb durch/Distribuito da:

And in the USA/Y en EEUU/En aux Etats-Unis/Und in USA/E in USA:

Denver, PA 17517 USA


Rev. 2014/04


0473

Sintering temperature:


1-5 units (not including pontics)

2-5 units (including pontics)

6-9 units:

Min°C: 1.450

Zirlux ST1

ZIRCONIUM DISCS / BLOCKS

INSTRUCTIONS FOR USE

1. Fix the blank / block according to the CAM-System instructions. Set the enlargement-factor into the unit..
2. Start the milling unit.
3. After milling, extract the blank / block with milled restoration.
4. Remove the restoration from the attached blank / block with a fissure bur.
5. Remove any dust from the restoration with pressured air before sintering.
6. Position the restoration on the sinter bed on the furnace tray.
7. Programme your furnace according to sintering instruction applicable for Zirlux ST1 shown in these instructions.
8. Run the sintering cycle.
9. Wait until cool down is completed before removing the sintered restoration.
10. Inspect the restoration in respect of flaws, wall thickness, seat and marginal fit. If necessary, make small adjustments with a water-cooled diamond bur.
11. Rinse the framework/anatomical crown with water and dry it. It is now ready for veneering. The framework is now ready for veneering /anatomical crown is ready for glazing.!

Safety precautions

- Design the restoration to maximize cross-sections.
- Carefully read the MSDS.

Sinter instructions

	Ramp Temp. °C / h	T. max °C	Dwell Time h	Cooling
Crowns	600	1450-1550	2	Natural furnace cooling
Bridges	300	1450-1550	2	

Indications for use

Zirconium discs / blocks are indicated for use as anatomy crown and a substructure for porcelain fused ceramic fixed dental restorations. Limitations are listed in table 1.

Table 1: Indications of use and maximum number of pontics.

Process chain	Material	Single unit crowns		Bridges, number of pontics	
		Anterior	Posterior	Anterior	Posterior
Porous blank/block is machined to enlarged framework/autonomic crown, then sintered and veneered	3Y-TZP	X	X	3	2

*X indicated; digits show maximum number of pontics. Limitation of Bridge, number of Pontics - All Zirlux-blank / blocks are solely for use by or on the order of a dental professional. They are not for use by the general public or over-the-counter.

Contraindications

Contraindications are all applications not indicated in table 1 above.

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