Calcium Hydroxide Baseline SC (GB)

Instructions for use

Product description:

DE Calcium hydroxide is a self-curing radiopaque calcium hydroxide paste for (in)direct pulp capping and for use as a liner under all normal filling materials.

Indications:

DE Calcium hydroxide is indicated for deep cavities to reduce the pulp reaction on sensitivity (isolation layer). In direct proximity of the pulp, the application of an aqueous calcium hydroxide preparation is recommended.

Application:

1. Mix equal quantities of base and catalyst (mixing ratio 1.2:1 g/g) on a mixing pad or glass plate for approx. 10 seconds.
2. The mixture should achieve a homogeneous consistency and uniform colour. Do not mix for too long!
3. Once mixed, insert the Calcium hydroxide material into the cavity as quickly as possible with a suitable instrument.

Working time:

Mixing time 10 seconds.

Working time 2 minutes (15 - 23°C).

Higher temperature than above mentioned and humidity accelerate the setting time and reduce the working time, causing the material to set more quickly in the patient's mouth.

Precaution:

DE Calcium hydroxide contains calcium hydroxide, salicylate, calcium wolframate, N-ethyltoluol-sulfonamide: do not apply in cases of sensitivities or allergies against any component.

Let totally cure before filling in order to reach a sufficient stability for the filling material.

Do not use the self curing calcium hydroxide as temporary filling or as the only liner carrying masticatory forces. The catalyst hardens at low temperatures.

The material has a strong alkaline effect, absolutely avoid contact with eyes (mucosa). Avoid contact with etchants (solubility).

Storage:

Store DE Calcium hydroxide at room temperature (15°C - 23°C). Reclose the tubes carefully after use. Be certain to match the colour-coded caps to the proper tube. Avoid humidity.

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DE Healthcare Ltd