# Nanohybrid Composite

## Nano-Hybrid Composite Resin-based dental restorative material **INSTRUCTIONS FOR USE**

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DE Healthcare Products  
Gillingham ME8 OSB U.K.

# NANOHYBRID COMPOSITE • NANO HYBRID COMPOSITE

## Resin based dental restorative material

### INSTRUCTIONS FOR USE

**Nanohybrid Nano Composite is a syringeable, sculptable, radiopaque, highly filled nano-hybrid composite for the restoration of all teeth. Nanohybrid Nano Composite covers a wide range of tooth shades with fewer shades. Nanohybrid Nano Composite, transparent is a highly translucent surface layer composite to improve gloss or to restore incisal edges. Nanohybrid Nano Composite bleach is useful for the restoration of whitened teeth.**

#### Composition Nanohybrid Nano Composite Duo Shade

• Methacrylates • Barium glass, silanized • Amorphous silica, hydrophobed

#### Technical Data

• Complies with ISO 4049:2000 • Average filler particle size: 0.6 µm • Range of particle size: 0.04–2.9 µm

#### Nanohybrid Nano Duo Shade

• Percentage by volume of total inorganic filler: 58 % • Percentage by weight of total inorganic filler: 77 %

#### Indications

Nanohybrid Nano Composite Duo Shade is specifically indicated for:

- Direct filling of all teeth
- More translucent material in the layer technique
- Bonding of composite and ceramic restorations

Nanohybrid Nano Composite, transparent is specifically indicated for:

- Translucent surface layers in all classes
- Reconstruction of incisal edges
- Shape and color corrections to improve esthetics
- Bonding of composite and ceramic restorations

Nanohybrid NanoComposite Bleach is specifically indicated for:

- Restoration of whitened teeth
- Restoration of deciduous teeth
- Characterization of chalk spots
- Masking of discolored teeth

#### Contra-indications

If the site cannot be isolated after enamel etching and during application and curing of Nanohybrid Nano Composite. If allergic to components of Nanohybrid Nano Composite. If oral hygiene is poor.

#### Side effects

No systemic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

#### Interactions with other agents

Agents which contain eugenol and/or oil of cloves may affect the polymerization of Nanohybrid Nano Composite. The use of zinc-oxide-eugenol cements in combination with Nanohybrid Nano Composite should therefore be avoided. Discolorations may occur when using cationic mouth rinses, plaque indicators or chlorhexidine.

#### Application - syringe

Cleansing the tooth. The tooth to be treated and the adjacent teeth are cleansed using a brush and prophylaxis-paste which does not contain fluoride. Selection of shade. Select shades before isolation. Ideally, the shade is determined in broad daylight using the VITA™ Lumin-Vacuum shade guide (registered trademark of the company VITA™). Shades are identical to those of the shade guide only 24 hours after curing. Isolation of the cavity (rubber dam) Adequate isolation is absolutely necessary for achieving optimum results. Preparation of the cavity: When preparing the cavity, every effort should be made to use a technique that is gentle to the tooth substance (principle of Adhesive Restoration Technique). Prepare the enamel and dentine using 80 µm preparation diamonds and finish using 25 µm finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material, thereby optimizing marginal seal. Pulp Protection: Covering the areas close to the pulp with hardsetting calcium hydroxide cements using the spot technique will provide protection of the pulp against bacterial infiltration. Bonding system: use a bonding according to the corresponding instructions for use.

Matrix and interdental wedges Apply a thin matrix band for class II, III and IV restorations. Fix the matrix proximally using interdental wedges. Application of Nanohybrid Nano Composite from a syringe. Dose Nanohybrid Nano Composite directly onto the mixing pad, take up in small increments with a composite instrument or any other application instrument, and apply to the cavity.

**Application – tip:** Remove the cap. Insert the composite capsules through the barrel opening of the dispenser. Rotate the angled composite capsule into the desired position. Extrude Nanohybrid Nano Composite into the cavity or onto an instrument by exerting gentle, even pressure on the handle. Working time Nanohybrid Nano Composite is photosensitive, and should not be exposed to light too long, in particular to the operating light or sunlight. If possible, protect exposed Nanohybrid Nano Composite material by means of light opaque covering (not blue).

#### Polymerization

Nanohybrid Nano Composite can be cured by using all current dental light sources.

Recommended exposure times:

<b>Curing time of layer Thickness: 2 mm</b>	$\leq 800 \text{ mW/cm}^2$	$> 800 \text{ mW/cm}^2$
Nanohybrid Duo Shade A1/D2	30 s	20 s
Nanohybrid Duo Shade A2/B2	30 s	20 s
Nanohybrid Duo Shade A3/D3	30 s	20 s
Nanohybrid Duo Shade A3.5/B3	40 s	20 s
Nanohybrid Duo Shade A4/M5	40 s	20 s
Nanohybrid Transparent	20 s	20 s
Nanohybrid Bleach	30 s	20 s

Nanohybrid Nano Composite must be cured layer by layer. Do not remove the inhibition layer since this would affect adhesion between the layers. In the event of underexposure there is a risk that the ideal strength values may not be achieved. Overexposure is not possible. Equipment whose light output is unknown should be checked with a light meter according to corresponding instructions for use.

### Finishing procedures

Finishing consists of three stages: coarse finishing, fine finishing and polishing. Finishing diamonds (40  $\mu\text{m}$  and 15  $\mu\text{m}$ ) are the most versatile and the least destructive. They should be used with very light pressure, with a constant wiping motion and a generous water spray. The ideal speed ranges between 5000 and 15 000 rpm. Flexible diamond files (80  $\mu\text{m}$ , 40  $\mu\text{m}$  and 15  $\mu\text{m}$ ) and strips are suitable for approximal surfaces. For high gloss polishing results use any comparable instruments. Fluoridize all treated surfaces after finishing.

### Emergency measures

Danger from swallowing of Nanohybrid Nano Composite can be excluded due to the small quantity applied and its lack of toxicity. In case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In the case of contact with the eyes, rinse thoroughly with water (10 min) and consult physician.

### Toxicological characteristics

In addition to the active components, Nanohybrid Nano Composite contains small amounts of stabilizers, activators and color pigments. Allergic reaction to resin components or the amine component of the activator system is possible.

### Note

Only supplied to dentists and dental laboratories or upon their instructions. Keep out of the reach of children! Do not mix different shades on the mixing pad, since this procedure will lead to air entrapment and porosities.

### Shelf life and marking

The expiry date and LOT number are marked on the syringes and tips.

### Storage

Store at 4–23 °C / 39–73 °F. Avoid exposure to direct sunlight or other heat sources.

### Caution

Federal law restricts this device to sale by or on the order of a dentist.

### Infection control

Compliance with the following instructions ensures optimal hygienic conditions when using Nanohybrid Nano Composite:

- Nanohybrid Nano Composite Capsules are for single use only.
- Recap composite syringes after use.

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