FOR PROFESSIONAL USE ONLY

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 **DE Healthcare Products**  
DEDICATED TO EXCELLENCE

**0473**

Distributed by:  
Distribuido por:  
Distribue par:  
Vertrieb durch:  
Distribuito da:

**DE Healthcare Products**  
Gillingham ME8 OSB U.K.

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## ZINC OXIDE EUGENOL TEMPORARY CEMENT

### Description:

**Zinc Oxide Eugenol Temporary Cement** is used for provisional cementation. Zinc Oxide Eugenol Temporary Cement contains eugenol with its associated properties. Do not use Zinc Oxide Eugenol Temporary Cement with resin bonding or place next to resin endodontic sealers.

### Automix Procedure:

1. Remove cap. Express a small amount on pad to ensure even flow.
2. Attach mixing tip.
3. Deliver to provisional.

**Zinc Oxide Eugenol Temporary Cement** set time: 3 - 3.5 minutes. Add a drop of water to accelerate set.

**Important Note:** *Polycarboxylate cement is acidic when initially placed. Patients should be anaesthetised or they could experience discomfort.*

### Manual Mix Procedure:

1. Remove cap. Express a small amount on pad to ensure even flow.
2. Mix conventionally with mixing pad and disposable spatula.
3. Deliver to provisional.

### Removal of Temporary Cement:

1. Remove provisional.
2. Scrape off residual cement with appropriate hand instrument.
3. Scour preparation.
4. Wash thoroughly.

### Precautions:

1. Do not use on patients with known sensitivity to polycarboxylate (polyacrylic acid) or to Eugenol.
2. While being stored or mixed, chemical cannot come in contact with metal otherwise polymerization process will begin.
3. Patients should be anaesthetised to prevent discomfort during cementation as polycarboxylate, and other acid based cements (e.g. glass ionomer and zinc phosphate) can cause discomfort at least until fully set.

Store at room temperature. **DO NOT REFRIGERATE.**

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