DE Amalgam SPHERICAL

High-copper, non gamma II
spherical amalgam alloy
INSTRUCTIONS FOR USE

CONTENTS:

Ag 60%, Sn 28%, Cu 12%

PROPORTIONS:

Powder: Alloy : mercury ratio 1:0.88

Capsules (Self Activated):

- 1 spill 400 mg alloy 290 mg mercury
- 2 spill 600 mg alloy 435 mg mercury
- 3 spill 800 mg alloy 580 mg mercury

Each self-activating capsule is ultrasonically welded for your protection.

TRITURATION:

Amalgamator	Speed	1 Spill	2 Spill	3 Spill
Silamat		4-5	5-7	7-9
Varimix	M2	4-5	5-7	7-9
Duomat (4000rpm)	3200	5-7	5-7	7-9
Capmix Espe		4-5	5-7	7-9
Wig-L-Bug	High (4600)	10-12	12-16	14-16

PLEASE DO NOT USE LOW SPEED AMALGAMATORS!

Not designed to be used with planetary mixers such as copernicus and capmix.

NOTE: The triturated amalgam should be bright and have a plastic consistency. If the amalgam is excessively plashy and wet looking and/or sets fast, decrease trituration time or amalgamator speed. A mix that appears excessively dry indicates that trituration time or amalgamator speed should be increased. Your clinical experience in obtaining a satisfactory working consistency will be the best guide to trituration times. Immediately after trituration – Firmly grasp bottom half of capsule with index finger and hold securely. Insert opening-rod into the indentation on the top of the capsule. Snap capsule open.

WORKING TIME:

Condense 2 - 3.5 min.

Band removal, carving and burnishing 3.5 - 7 min.

Polish 24 hr. after placement.

SPECIAL NOTE

1. This alloy requires less mercury than other alloys, when triturating with Dentomat adjust the regulating screw for a leaner mixture.

2. Refer to American Dental Association Publication, "Recommendations in Dental Mercury Hygiene"

TEST CONDITIONS

Testing to A.D.A. Specification No. 1 and B.S. 21559:1992

Alloy/Mercury ratio : 1:0.88

Weight of alloy : 0.6 gr

Amalgamator time (for ordinary units) : 4-5 sec.

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