TEMPORARY CROWN & BRIDGE MATERIAL

Material provisional para coronas y puentes
Matériau pour couronnes et bridges provisoires
Provisorisches Kronen- und Brückenmaterial
Materiale per corone e ponti provisorio
Noodkroonmateriaal

9792680 A2 starter kit
9792681 A3 starter kit

DEHP™ DE Healthcare Products
DEDICATED TO EXCELLENCE

UK Instructions for Use

Temporary Crown & Bridge material is a fluorescent cold polymerising paste-paste system for the production of temporary crowns, bridges, inlays and onlays. **Temporary Crown & Bridge material** consists of base paste and catalyst paste.

Application

Temporary Crown & Bridge material in the cartridge:

Insert **Temporary Crown & Bridge material** cartridge into a suitable dispensing system. There may be a slight difference in the levels in each chamber. Remove the cap and dispense the material until it is evenly distributed from each chamber. Attach the mixing tips, ensuring that the notches are lined up and turn clockwise until the tip locks. Extrude the material through the mixing tip. After use, do not remove the mixing tip, this will act as the cap. Change the tip when you are ready to use the product again. Check that the material can flow correcty from the cartridge.

Directions for Use:

Level out any undercuts in the preparation using a non-eugenol cement.

Prepare an impression if not using prefabricated crowns or similar restoratives. Prepare the impression by levelling undercuts and cutting drainage channels.

If the wall thickness of the temporary restoration is too thin, enlarge the impression, in particular remove the septums in the interdental embrasures. Then clean and air-dry the impression.

Apply **Temporary Crown & Bridge material**directly into the impression. Do not overfill.

Temporary Crown & Bridge material reaches an elastic phase after 1 - 1.5 min from beginning of mixing (approx. 0.5 - 1 min being in the mouth).

The temporary should be removed from the mouth during this elastic phase (with 1.5min of mixing)

The setting process can also be monitored by observing the excess material inthe mouth.

The setting time of residual material n the mixing tips cannot be compared with the setting time of the material in the mouth.

At a room temperature of 23°C and 50% humidity working times are as follows:

min:s 0:00	min:s 0:30	min:s 1:00	min:s 1:30	min:s 4:00	
application in the mouth	setting in the mouth	removal setting	complete	finishing	

Slight differences may result due to different room temperatures.

After approx. 2 - 3 min take the temporary out of the impression and wipe dry (if necessary, use a solvent, e.g. ethyl alcohol). This Inhibition [smear] layer usually caused by oxygen and allows easier finishing. The surface of the temporary restoration should appear silky-matt.

Carefully trim any excess material and roughen areas near the gingiva. Cement the temporary restoration using eugenol-free temporary cement.

Indications/precautions:


- Temporary Crown & Bridge material**contains methacrylates, amines, terpenes, benzoylperoxide, BHT. **Temporary Crown & Bridge material** should not be used in case of known hypersensitivities (allergies) to these ingredients.
- Before making the temporary, the prepared cores can be insulated for desensitisation (fluoride varnishes, core varnishes, etc.).
- In order to avoid fractures or injuries from debris the temporary should not be subjected to excessive loads (molar area).
- Linings and cores built up with resin-based materials (composites) should be sufficiently insulated (e.g. glycerine) to avoid removal with the temporary.
- Inadequate hygiene or intensive use of certain mouth rinses may lead to discolouration.
- Finishing can be completed once the material is fully set (approx. 4 mins after mixing).
- Restoration of defective temporaries:
- In the RARE occurrence that the temporary breaks the following procedures **Temporary Crown & Bridge material**, the following procedures are recommended:
 - Air bubbles: fill with **Temporary Crown & Bridge material**
 - In case of fracture immediately after the manufacturing of the temporary: re-connect the fracture with freshly mixed **Temporary Crown & Bridge material**
 - Restoration of already used temporaries: Roughen the fracture with a bur and place undercuts. Apply bond on the prepared fractures and cure according to the respective instructions for use. Apply a flowable composite to the prepared areas. Reassemble the fractured pieces and cure according to the instructions for use.

Storage

Store between 4°C to 23°C. Do not use after expiry date.

Protect from bright light, i.e. store in closed packages, drawers, etc.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/ or written information is given to the best of our knowledge and without obligation. Our information and/or advice does not relieve you from examining the materials for the suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

Additional products: <i>9790316</i> Universal Automix Dispenser 1:1 <i>9792409</i> Mixing canulas 1:1, light blue, 50 pcs										
9792680	A2 Starter kit	1		6		blue 1:1				
9792681	A3 Starter kit	1		6		blue 1:1				

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