DE Healthcare Ltd
Unit 9, Kingsthorpe Business Centre,
Studland Road, Kingsthorpe,
Northampton NN2 6NE U.K.

Distributed by/Distribuido por/Distribué par/Vertrieb durch/Distribuito da:

And in the USA,/Y en EEUU,/Et aux Etats-Unis,/Und in USA,/E in USA:

41 Weaver Road
Denver, PA 17517 USA

Rev. 2014/05

SINGLE BOND 7TH GENERATION

Universal adhesive • self-etching, single component light-cured

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
DE Healthcare Products
DEDICATED TO EXCELLENCE

UK SINGLE BOND 7TH GENERATION

IS A SINGLE COMPONENT LIGHT-CURED, SELF-ETCHING ADHESIVE FOR ADHESIVE RESTORATION TECHNIQUES.

INSTRUCTIONS FOR USE

Composition

Single Bond 7th Generation contains:

Methacrylates
Photo initiators
Ethanol
Water

Indications for use

Adhesive bonding for use with adhesive restorative techniques:

- Adhesive bond of composite (e.g. Nanocomposite from DEHP*) and compomer materials to enamel and dentin
- Adhesive bond of ceramic and composite restorations (e.g. Nanocomposite from DEHP*) to enamel and dentin
- Adhesive bond of composite material to ceramics, composite material (e.g. Nanocomposite from DEHP*), metal and amalgam (e.g. Dental Amalgam from DEHP*)
- Dentin sealing

Contraindications

Hypersensitivity due to any of the elements within DEHP* Single Bond 7th Generation. Inadequate oral hygiene. If the working area cannot be kept completely dry during application.

Side effects

Elements of Single Bond 7th Generation may lead to sensitivity or an allergic reaction for patients with a predisposed condition.

Interactions

Phenolics and other substances (e.g. zinc oxide eugenol) may not come in contact with Single Bond 7th Generation, as they will inhibit polymerisation.

Application**Note**

To prevent the bond from curing in the dispensing well, Single Bond 7th Generation should be immediately covered with the black cover after each application.

1. Conventional filling technique**1.1 Cavity preparation**

Cavity preparations should be prepared according to the principles of the adhesive restoration technique. If indicated, apply hard-setting calcium hydroxide cement (e.g. DEHP* Calcium Hydroxide) in areas near the pulp, using the spot over-capping technique. Prepare the enamel using a finishing bur. In cases with enamel that is not prepared, etch it using phosphoric acid (e.g. DEHP* Etch Gel) or cleanse it using pumice. Thoroughly rinse and dry the enamel afterwards. For hygienic reasons, it is recommended to use a dental dam (e.g. 9792268 DEHP* Latex Dental Dam or 9792146 DEHP* Non-Latex Dental Dam).

1.2 Application

Shake the bottle well before use and dispense one drop into the dispensing well. Massage Single Bond 7th Generation using a disposable brush for 20 s onto the dentin and enamel. Gently dry for 5 s using oil-free compressed air.

Light-cure for 10 s (e.g. Halogen- or LED Polymerisation units with a light output of $>800 \text{ mW/cm}^2$). If using a curing lamp with a light output less than this, be sure to increase the length of the curing time accordingly.

1.3 Composite and compomer fillings

Apply restorative material according to the manufacturer`s instructions for use.

2. Adhesive bonding with inlays/onlays, veneers and crowns

The application of Single Bond 7th Generation is only possible under the condition that sufficient light access (according to the manufacturer`s instructions for use) can be achieved through the complete restoration.

2.1 Preparation of the restoration

The preparation of the restoration is dependent on the material used and according to the manufacturer`s instructions for use.

2.2 Cavity Preparation

Proceed as described in point 1.1.

2.3 Application

Proceed as described in point 1.2. Avoid over-application of the adhesive as this could impact the fit of the restoration.

2.4 Cementing

According to the manufacturer`s instructions for use.

3. Adhesive bonding with other surfaces

3.1 Preparation

3.1.1 Ceramics

Intraoral application:

Roughen ceramic surface using a diamond bur or sandblasting, etch with 35 % phosphoric acid, rinse, and dry using oil-free compressed air.

Extraoral application:

Etch ceramic surface using 5 % HF-acid, rinse, and dry using oil-free compressed air.

Optional:

Ceramic surfaces can also be silanised (according to the manufacturer's instructions for use).

3.1.2 Metal and amalgam

Roughen metal surfaces using a diamond bur or sandblasting. Clean using steam jet (extraoral) or acetone; and dry using oil-free compressed air.

3.1.3 Composite

Roughen composite surfaces using a diamond bur or sandblasting, etch with phosphoric acid, rinse and dry using oil-free compressed air.

3.2 Application

Shake the bottle well before use and dispense one drop into the dispensing well. Afterwards close the bottle. Massage Single Bond 7th Generation using a disposable brush for 20 s onto the surface. Gently dry for 5 s using oil-free compressed air.

Light-cure for 10 s (e.g. Halogen- or LED Polymerisation units with a light output of >800 mW/cm²). If using a curing lamp with a light output less than this, be sure to increase the length of the curing time accordingly.

3.3 Composite and compomer fillings

Apply restorative material according to the manufacturer's instructions for use.

4. Dentin sealing

Proceed as described in points 1.1 and 1.2.

Emergency measures

In case of direct contact with oral mucosa, it is sufficient to rinse area using tap water. In case of contact with eyes, rinse eyes thoroughly using tap water. Seek immediate professional medical attention from an ophthalmologist or eye care specialist.

Notes

For use only by dentists and dental technicians, or upon their request. Keep out of the reach of children!
Only recycle completely empty containers.

Shelf life and marking

Expiry date and number are marked on the immediate containers and packaging.

Do not use after the expiry date.

Storage at temperature 2–8 °C / 35,6–46 °F.

Keep out of direct sunlight or near sources with warm temperatures.

Infection control

Compliance with the following recommendations will help ensure optimal hygienic conditions, when using Single Bond 7th Generation:

- Only use disposable brush for one patient.
- After each use, close the bottle with the cap, and disinfect using conventional methods.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

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Composicion

Single Bond 7th Generation se compone de:

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