# NANOHYBRID COMPOSITE FLOW

## INSTRUCTION FOR USE

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ISO 4049, Type 1, class 2, group 1, radio-opaque

DE Healthcare Ltd  
Unit 9, Kingsthorpe Business Centre,  
Studland Road, Kingsthorpe,  
Northampton NN2 6NE

Distributed by/Distribuido por/Distribuit per/Vertriebt durch/Distribuito da:  
And in the USA/En un EEUU/Et aux Etats-Unis/Und in USA/E in USA:

Denver, PA 17517 USA

Rev. 2013/11

# GB - NANOHYBRID COMPOSITE FLOW

## INSTRUCTIONS FOR USE

Nanohybrid Composite Flow is a nanofilled, radiopaque flowable composite for special indications. Nanohybrid Composite Flow is cured in a tooth cavity or on a model with all current light-curing lamps (e.g. BA light curing units).

### Composition

- Methacrylates Barium glass, silanized Amorphous silica, hydrophobic

### Technical Data

- Complies with ISO 4049:2000 • Average filler particle size: 0.6  $\mu\text{m}$  • Range of particle size: 0.04–2.5  $\mu\text{m}$

### Nanohybrid Composite Flow

- Percentage by volume of total inorganic filler: 42 % • Percentage by weight of total inorganic filler: 63 %

### Indications

- Direct fillings of class V (cervical caries, root erosion, wedge-shaped defects)
- Repairs to anterior teeth (class III, IV)
- Small fillings in the posterior region
- Fillings of minimal cavities
- Repairs of composite and ceramic veneers
- Block out of undercuts
- Adhesive cementing of composite and ceramic restorations, if light penetration is possible.
- Preventive resin restorations

### Contra-indications

If any allergies exist to any components within Nanohybrid Composite Flow. If the application field cannot be isolated during application and curing of Nanohybrid Composite Flow. For patients with poor oral hygiene.

### Side effects

No systematic side effects are known. Contact allergies with products of similar composition have been reported in isolated cases.

### Interactions with other agents

Agents containing eugenol and/or oil of cloves may affect the polymerization of Nanohybrid Composite Flow. The use of zinc-oxide eugenol cements in combination with Nanohybrid Composite Flow should therefore be avoided. Discolorations may occur when using cationic mouth rinses as well as plaque indicators or chlorhexidine.

**Note:**

- Light can lead to premature curing. Recap syringe after each use.
- To prevent run off, pull back the syringe plunger.
- For repairs, the treated surface must be roughened first.

**General information****Preliminary preparation:****- Cleansing the tooth**

Cleanse tooth to be treated and the adjacent teeth using a brush and prophylaxis paste which does not contain fluoride.

**- Selection of shade**

Selection of shade should be made before isolation. The shade is ideally determined in broad daylight using the Nanohybrid Composite Shade Guide or the Vita™ Shade Guide. Shades reach their final tone 24 hours after curing.

**- Isolation of the cavity (rubber dam)**

Adequate isolation is absolutely necessary for achieving optimum results. The use of DE Healthcare Dental Dams is recommended.

**- Matrix and inter-dental wedges**

When working proximally, apply a thin matrix band. Fix the matrix proximally using interdental wedges (eg DE Healthcare wooden wedges).

**- Adhesive system**

Apply adhesive system (e.g. DE Healthcare bond) according to the corresponding instructions for use.

**- Working time**

Nanohybrid Composite Flow is light sensitive and should not be left too long under intensive light before curing, especially operatory light or sun light.

**- Polymerization**

Nanohybrid Composite Flow can be cured using all current curing lights (e.g. BA Light curing units).

**Recommended exposure times:**

<b>Layer Thickness: 2 mm</b>	> 500 mW/cm <sup>2</sup>
Nanohybrid Composite Flow A2/B2	20 s
Nanohybrid Composite Flow A3/D3	20 s

## **Finishing**

For the tooth preparation, appropriate rotary instruments and polishers. Fluoridate all surfaces after completion.

## **Application of Nanohybrid Composite Flow for direct restorations**

Preparation of the cavity

When preparing the cavity, every effort possible should be made to use a technique that preserves tooth structure (principle of the Adhesive Restoration Technique). Prepare enamel and dentin using 80  $\mu\text{m}$  preparation diamonds and finish using 25  $\mu\text{m}$  finishing diamonds. Bevelled enamel margins are recommended to increase the area of adhesion between tooth and filling material and there by optimize the marginal seal.

## **Pulp Protection**

Covering the areas close to the pulp with hard-setting calcium hydroxide cements using the spot technique protects the pulp against bacterial infiltration.

## **Application of Nanohybrid Composite Flow**

Apply the material using the mounted applicator tip directly into the cavity.

## **Application of Nanohybrid Composite Flow for indirect restorations**

Pretreat contact surfaces of cavity according to the corresponding manufacturer's instructions for use.

## **Application of Nanohybrid Composite Flow**

Dispense the material using the mounted applicator tip directly into the restoration or cavity preparation. Position restoration in place using gentle pressure.

## **Removal of excess material**

Roughly remove excess material (e.g. with a brush or spatula). Afterwards, hold restoration with pressure in position and remove additional excess material with a spatula or a suitable instrument.

## **Polymerization**

Recommended curing time: cure each side for 40 s.

## **Emergency measures**

An acute danger of Nanohybrid Composite Flow can be excluded due to the small quantity applied and the lack of toxicity. In the case of direct contact with the oral mucosa, rinsing with tap water is sufficient. In case of contact with the eyes, rinse thoroughly with water (10 min), and consult an ophthalmologist showing these instructions for use.

**Note**

Only supplied to dentists and dental laboratories or upon their instructions. Keep out of the reach of children! Should not be used after expiration date.

**Shelf life and marking**

The expiration date and the lot number are marked on the packaging containers.

**Storage**

4–23 °C / 39–73 °F. Avoid exposure to direct sunlight or other heat sources.

**Caution**

Federal law restricts this device to sale by or on the order of a dentist.

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