# Zirlux<sup>®</sup> Wax

Wax Discs

Disco in cera

Disque de cire

Wachs-Scheiben

Cera disco

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Ø 98,5 mm

H: 20 mm

DE Healthcare Ltd  
Unit 9, Kingsthorpe Business Centre,  
Studland Road, Kingsthorpe,  
Northampton NN2 6NE

Distributed by/Distribuido por/Distribué par/Vertriebt durch/Distribuito da:  
And in the USA/Y en EEUU/Et aux Etats-Unis/Und in USA/E in USA:

Denver, PA 17517 USA

Rev. 2014/01

**9793164**

Blue

Azul

Blue

Blau

Blu

Blauw

**9793165**

Grey

Gris

Gris

Gray

Grigio

Grijs

**9793166**

White


Blanco

Blanc

Weiß

Bianco

Wit

## INSTRUCTIONS for dentists/dental technicians

Use: For the manufacturing of the wax scaffolding of dental implants with plastic covered crowns, through CAD/CAM systems.

### 1. General

The wax discs can be processed with different milling systems.

### 2. Storage

Protect from humidity and sunlight.

### 3. Area of application:

The forming wax is not designed for use in oral application. The material can be used for the preparation of dental crowns and bridges. Since it burns out leaving without residues, it is suitable for the preparation of casting moulds and ceramic moulding processes.

### 4. Processing

The wax milled scaffoldings can be removed with the help of appropriate cutters or cutting discs.

### 5a. Installation of the casting channel

The setting and the selection of the diameter of the casting channel depends on the used devices and alloys or pressed ceramics.

### 5b. Embedding

Line the silicon or metal muffle with a liner. Mix and embed the investment material under vacuum observing the instructions provided by the manufacturers of the investment material and alloys.

### 6. Disposal

Uncontaminated residues can be re-used for wax related activities. Since the product is insoluble in water, it can be disposed of together with normal household waste.

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